

Polyclonal Guinea-Pig anti-Glucagon

Clone no. -

MONOSAN

Product name	Polyclonal Guinea-Pig anti-Glucagon
Host	Guinea pig
Applications	IHC-fr, IHC-P (1:1,000-1:2,000)
Species reactivity	Cat, Human
Conjugate	-
Immunogen	Highly purified porcine pancreatic glucagon (Novo)
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	100 µl
Concentration	n/a
Format	-
Storage buffer	Contains 0.09% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Guinea-Pig anti-Glucagon

Clone no. -

MONOSAN

Additional info

Glucagon is a common constituent of endocrine pancreatic tumors and of rectal carcinoids. The antibody is specific to pancreatic glucagon.

Absorption with 10-100 ug glucagon per ml diluted antiserum abolishes the staining.

Positive control: Bouin-fixed paraffin sections of cat pancreas.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES