

Polyclonal Rabbit anti-Ornithine decarboxylase

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Ornithine decarboxylase
Host	Rabbit
Applications	IHC-fr (1:300-1:500)
Species reactivity	Human, Mouse, Rat
Conjugate	-
Immunogen	Recombinant human ornithine decarboxylase produced in E. coli
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Ornithine decarboxylase

Clone no. -

MONOSAN

Additional info

Ornithine decarboxylase is the rate-limiting enzyme in polyamine biosynthesis converting ornithine into putrescine. In normal tissue ornithine decarboxylase activity is low but increases in proliferating tissue.

Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes staining.

Positive control: Stefanini-fixed frozen sections of renal cortex from testosterone-treated mice.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES