

Product datasheet PS261

MONOSAN[®]

Polyclonal Rabbit anti-CEA

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-CEA
Host	Rabbit
Applications	IHC-fr, IHC-P (1:1,000-1:5,000)
Species reactivity	Human
Conjugate	-
Immunogen	Human carcinoembryonic antigen (CEA; Rougier Bio-Tech code 70015)
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-CEA

Clone no. -

MONOSAN

Additional info

The antibody reacts with a large variety of tumors with oncofetal characteristics; the antigen may also be detected in normal epithelia and tissue of non-neoplastic state.
Absorption with 10-100 ug immunogen per ml diluted antiserum abolishes staining.

Positive control: Bouin-fixed paraffin sections of human colon carcinoma.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES