Product datasheet

PS266



MONOSAN

Polyclonal Rabbit anti-VAChT Clone no. -

Product name	Polyclonal Rabbit anti-VAChT
Host	Rabbit
Applications	IHC-fr (1:800-1:1,200)
Species reactivity	Human
Conjugate	-
Immunogen	Synthetic peptide from the C-terminus of human VAChT (11 amino acid residues) conjugated to BSA
lsotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS266



Polyclonal Rabbit anti-VAChT Clone no. -

MONOSAN

Additional info

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com