Product datasheet	PS279
-------------------	-------

MONOSAN

Polyclonal Rabbit anti-Csk-binding protein (PAG) Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Csk-binding protein (PAG)
Host	Rabbit
Applications	IP, WB, IHC-P
Species reactivity	Mouse, Human
Conjugate	-
Immunogen	Recombinant intracellular fragment (aa 97-432) of human Cbp (PAG).
lsotype	-
Clonality	Polyclonal
Clone number	-
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS279

Polyclonal Rabbit anti-Csk-binding protein (PAG) Clone no. -

MONOSAN[°]

MONOSAN

Additional info

PAG (phosphoprotein associated with GEMs), also known as Cbp (Csk-binding protein), is a ubiquitously expressed 46 kDa transmembrane adaptor protein present in membrane rafts (glycosphingolipid-enriched microdomains), which however migrates on SDS PAGE gels anomalously as an 80 kDa molecule. Following tyrosine phosphorylation by Src family kinases, PAG binds and thereby activates the protein tyrosine kinase Csk, the major negative regulator of the Src family kinases. Signaling via the B-cell receptor in B cells or high affinity IgE receptor (FcepsilonRI) in mast cells leads to PAG increased tyrosine phosphorylation, loss of Csk binding and increased activation of the protein tyrosine kinase Lck.

References	1.	-
	2	-
	3.	-
	4.	-
	5.	-

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com