

Product datasheet PS281F

MONOSAN[®]

Polyclonal Sheep anti-C3c, FITC

Clone no. -

MONOSAN

Product name	Polyclonal Sheep anti-C3c, FITC
Host	Sheep
Applications	IHC-fr
Species reactivity	human
Conjugate	FITC
Immunogen	Human C3 purified from serum.
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	10 mg/ml
Format	-
Storage buffer	PBS
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Sheep anti-C3c, FITC

Clone no. -

MONOSAN

Additional info

Sheep anti Human C3c antibody recognizes the C3c component of human complement, formed as a result of the inactivation of C3b. Sheep anti Human C3c antibody may be used for the detection of C3 deposits in tissues following complement activation.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES