

Product datasheet PS282F

MONOSAN[®]

Polyclonal Sheep anti-C1q, FITC

Clone no. -

MONOSAN

| | |
|---------------------------|----------------------------------|
| Product name | Polyclonal Sheep anti-C1q, FITC |
| Host | Sheep |
| Applications | IHC-fr |
| Species reactivity | human |
| Conjugate | FITC |
| Immunogen | Human C1q, purified from plasma. |
| Isotype | - |
| Clonality | Polyclonal |
| Clone number | - |
| Size | 1 ml |
| Concentration | 10 mg/ml |
| Format | - |
| Storage buffer | PBS |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS282F

MONOSAN[®]

Polyclonal Sheep anti-C1q, FITC

Clone no. -

MONOSAN

Additional info

Sheep anti Human C1q antibody recognizes human C1q, a sub-component of the C1 complement component.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES