Product datasheet PS290



Polyclonal Rabbit anti-Mycobacterium tuberculosis

Clone no. - MONOSAN

Product name Polyclonal Rabbit anti-Mycobacterium tuberculosis

Host Rabbit

Applications IHC-P,IHC-fr

Species reactivity bacterial

Conjugate Purified

Immunogen Purified protein derivative (PPD) from <i>M. tuberculosis</i>

Isotype -

Clonality Polyclonal

Clone number -

Size 1 ml

Concentration 4.0 mg/ml

Format -

Storage buffer PBS with 0.1% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS290



Polyclonal Rabbit anti-Mycobacterium tuberculosis

Clone no. - MONOSAN

Additional info

Rabbit anti Mycobacterium tuberculosis polyclonal antibody recognizes PPD from Mycobacterium tuberculosis. Rabbit anti M. tuberculosis has not been cross absorbed and may react with related micro-organisms, however the antibody is non-reactive with E.coli K12, Salmonella typhimurium, Pseudomonas aeruginosa, Streptococcus (group B), Candida albicans and Neisseria meningitidis.

References 1.

2 -3. -

4. -5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES