Product datasheet PS308

MONOSAN

Rabbit anti-CD3

Clone no. - MONOSAN

Product name Rabbit anti-CD3

Host Rabbit

Applications IHC-P (1:100-1:1000)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

Clonality Polyclonal

Clone number -

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS308

MONOSAN

Rabbit anti-CD3 Clone no.

MONOSAN

Additional info

Anti-CD3 antibody has been considered the best all around T-cell marker. This antibody reacts with an antigen present in early thymocytes. The positive staining of this marker may represent a sign of early commitment to the T-cell lineage.

References 1.	Beverley PC, et al. Eur	J Immunol. 1981; 11:3
---------------	-------------------------	-----------------------

2 Clevers H, et al. Eur J Immunol. 1988; 18:705-10

3. Hedvat CV, et al. Hum Pathol. 2002; 33:968-74

4. Karube K, et al. Am J Surg Pathol. 2003; 27:1366-74

5. Dogan A, et al. Am J Surg Pathol. 2003; 27:903-11

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES