Product datasheet PS309

**MONOSAN** 

Rabbit anti-CEA

Clone no. - MONOSAN

Product name Rabbit anti-CEA

**Host** Rabbit

Applications IHC-P (1:200-1:500)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

**Clonality** Polyclonal

Clone number -

Size 1 ml

**Concentration** n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

# Product datasheet

PS309

**MONOSAN** 

Rabbit anti-CEA Clone no.

MONOSAN

### Additional info

Anti-CEA specifies a group of proteins in the Carcinoembryonic Antigen (CEA) family of proteins which are present in the epithelia of various types and tumors (both benign and malignant) derived from such epithelia. Such tissues are represented by the epithelia of colon, bronchus, alveoli, breast, pancreas, biliary tract, superficial layer and parietal layers of the stomach. Predominately biliary canaliculi are labelled in the liver and this factor is useful in the diagnosis of hepatocelluar carcinoma. Anti-CEA has been quite useful in differentiating adenocarcinoma of the lung vs. mesothelioma. Associated products: CK 5/6, Calretinin, WT-1, E-Cadherin, TTF-1, TAG-72, EMA, CK 20

# References

- 1. Shield PW, et al. Am J Clin Pathol. 1996; 105:157-62
- 2 Sheahan K, et al. Am J Clin Pathol. 1990; 94:157-64
- 3. Alkushi A, et al. Virchows Arch. 2003; 442:271-7
- 4. Morrison C, et al. Mod Pathol. 2002; 15:1279-87
- 5. -

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES