Product datasheet

P5313	
1 2212	



MONOSAN

Polyclonal antibody against FSH Clone no. -

Product name	Polyclonal antibody against FSH
Host	Rabbit
Applications	IHC-P (1:25-1:100)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS313



Polyclonal antibody against FSH Clone no. -

Additional info

Anti-FSH is a useful marker in classification of pituitary tumors and the study of pituitary disease. It reacts with FSH-producing cells (gonadotrophs).

	Refe	rences	
--	------	--------	--

1.

Baenziger JU, et al. Biochim Biophys Acta. 1988; 947:287-306

- 2 Nussey SS, et al. BIOS Scientific Publishers Ltd; 2001 p. 217-79
- 3. Uccella S, et al. Pituitary. 2000; 3:131-9
- 4. Schmid M, et al. Pathol Res Pract. 2001; 197:663-9
- 5. La Rosa S, et al. Virchows Arch. 2000; 437:264-9

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com