## Product datasheet



MONOSAN

Polyclonal antibody against Gastrin Clone no. -

Product name	Polyclonal antibody against Gastrin
Host	Rabbit
Applications	IHC-P (1:100-1:500)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Product datasheet

Polyclonal antibody against Gastrin Clone no. -

# **MONOSAN**<sup>°</sup>

MONOSAN

## Additional info

Anti-Gastrin antibody gives positive staining of G-cells of human antral/pyloric mucosa and cells producing gastrin or a structural gastrin analogue as is seen in stomach; no staining of other cells or tissue types has been observed. This antibody may react with sulfated and non-sulfated forms of gastrin. The antibody cross-reacts with more than 50% of the present choleocystokinin octapeptide.

References	1.	Kasacka W, et al. Folia Morphol. 2012; 71:39-44.
	2	Hur K, et al. J Cancer Res Clin Oncol. 2006; 132:85-91
	3.	Waldum et al. Frontiers in Endocrinology. 2017; 8:1-7
	4.	-
	5.	-

# FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com

FOR-044 22-12-2021