Product datasheet PS315



Polyclonal antibody against GH

Clone no. - MONOSAN

Product name Polyclonal antibody against GH

Host Rabbit

Applications IHC-P (1:100-1:500)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

Clonality Polyclonal

Clone number -

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS315



Polyclonal antibody against GH

Clone no. - MONOSAN

Additional info

Anti-GH is a useful marker in classification of pituitary tumors and the study of pituitary disease (acromegaly). It reacts with GH-producing cells. Growth hormone receptors have been found in various non-pituitary cells, including that from hepatocellular carcinoma and various benign and malignant cutaneous lesions.

References 1. Rezaei M, et al. J Res Med Sci. 2012; 17:681-5

2 Al-Brahim NY, et al. J Clin Pathol. 2006; 59:1245-53

3. Fukaya T, et al. Cancer. 1980; 45:1598-1603

4. Kovacs K, et al. Virch Arch Pathol Anat. 1982; 395:59-68

5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES