Product datasheet PS319



Polyclonal antibody against hCG

Clone no. - MONOSAN

Product name Polyclonal antibody against hCG

Host Rabbit

Applications IHC-P (1:50-1:250)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

Clonality Polyclonal

Clone number -

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS319



Polyclonal antibody against hCG

Clone no. - MONOSAN

Additional info

hCG is a protein secreted in large quantities by normal trophoblasts; the antibody detects cells and tumors of trophoblastic origin such as Choriocarcinoma. Large Cell Carcinoma and Adenocarcinoma of Lung demonstrate hCG positivity in 90% and 60% of cases respectively. 20% of Squamous Cell Lung Carcinomas are positive for hCG. hCG expression by nontrophoblastic tumors may indicate aggressive behavior since it has been observed that hCG may play a role in the host response to a given tumor.

References 1.

- 1. Hes O, et al. Ann Diagn Pathol. 2014; 18:89-94
- 2 Mazina O, et al. Sci. Rep. 2017; 7:42219
- 3. -
- 4.
- 5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES