Product datasheet PS327



Polyclonal antibody against IgG

Clone no. - MONOSAN

Product name Polyclonal antibody against IgG

Host Rabbit

Applications IHC-P (1:100-1:1500)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

Clonality Polyclonal

Clone number -

Size 1 ml

Concentration n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS327



Polyclonal antibody against IgG

Clone no. - MONOSAN

Additional info

Anti-lgG reacts with immunoglobulin gamma (lgG) chains. This antibody is useful when differentiating and sub-classifying hematolymphoid neoplasms.

References	1.	Arnold A,	et al. New Eng J	Med. 1983	: 309:1593-	1599
------------	----	-----------	------------------	-----------	-------------	------

- 2 Leong AS, et al. Geenwich Medical Media Ltd. 1999; 217-219
- 3. Ando K, et al. Intern Med. 2000; 39:170-5
- 4. Taylor CR. Arch Path Lab Med. 1978; 102:113-121
- 5. Warnake R, et al. Masson Publishing USA. 1981; 203-221

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES