

Product datasheet PS330

MONOSAN[®]

Polyclonal antibody against LH

Clone no. -

MONOSAN

Product name	Polyclonal antibody against LH
Host	Rabbit
Applications	IHC-P (1:100-1:500)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against LH

Clone no. -

MONOSAN

Additional info

Luteinizing hormone (LH) is a heterodimeric glycoprotein produced by gonadotropic cells of the pituitary gland. Anti-LH is a useful marker to aid in the classification of pituitary tumors and the study of pituitary disease.

References

1. Sano T, et al. Virchows Arch A Pathol Anat Histopathol. 1990; 417:361-7
2. Felix I, et al. Hum Pathol. 1991; 22:719-21
3. Saccomanno K, et al. J Clin Endocrinol Metab. 1994; 78:1103-7
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES