

Product datasheet PS334

MONOSAN[®]

Polyclonal antibody against Myoglobin

Clone no. -

MONOSAN

Product name	Polyclonal antibody against Myoglobin
Host	Rabbit
Applications	IHC-P (1:50-1:200)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against Myoglobin

Clone no. -

MONOSAN

Additional info

Immunostaining with anti-myoglobin provides a specific, sensitive, and practical procedure for the identification of tumors of muscle origin. Since myoglobin is found exclusively in skeletal and cardiac muscle and is not present in any other cells of the human body, it may be used to distinguish rhabdomyosarcoma from other soft tissue tumors.

References

1. Mukai K, et al. Am J Surg Pathol. 1979; 3:373-6
2. Corson JM, et al. Am J Pathol. 1981; 103:384-9
3. Brooks JJ. Cancer. 1982; 50:1757-63
4. Furlong MA, et al. Ann Diagn Pathol. 2001; 5:199-206
5. Kahn HJ, et al. Cancer. 1983; 51:1897-903

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES