Product datasheet PS335



## Polyclonal antibody against Napsin A

Clone no. - MONOSAN

Product name Polyclonal antibody against Napsin A

**Host** Rabbit

Applications IHC-P (1:100-1:500)

Species reactivity human

Conjugate -

Immunogen Unknown or proprietery to MONOSAN and/or its suppliers

Isotype -

**Clonality** Polyclonal

Clone number -

Size 1 ml

**Concentration** n/a

Format -

Storage buffer Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide

Storage until expiry date 2-8°C

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet PS335



Polyclonal antibody against Napsin A

Clone no. - MONOSAN

## Additional info

Napsin is a pepsin-like aspartic proteinase in the A1 clan of the AA clade of proteinases. There are two closely related napsins, napsin A (NAPSA) and napsin B (NAPSB). Napsin A is involved in processing propeptide pulmonary surfactant protein B (proSP-B) in the lung.4 In normal tissue, Napsin A is expressed in type II pneumocytes of the lung and proximal tubules of the kidney. Napsin A is a useful marker for lung adenocarcinoma

References 1. Brasch F, et al.   Biol Chem. 2003; 278: 490
--

- 2 Jagirdar J et al. Arch Pathol Lab Med. 2008; 132:384-96
- 3. Bishop JA, et al. Hum Pathol. 2010; 41:20-5
- 4. Ye J, et al. Appl Immunohistochem Mol Morphol. 2011; 19:313-17
- 5. Mukhopadhyay S, et al. Am J Surg Pathol. 2011; 35:15-25

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES