Product datasheet

PS341



MONOSAN

Polyclonal antibody against TdT Clone no. -

Product name	Polyclonal antibody against TdT
Host	Rabbit
Applications	IHC-P (1:100-1:500)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietery to MONOSAN and/or its suppliers
lsotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	_
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

PS341

MONOSAN

Polyclonal antibody against TdT Clone no. -

MONOSAN

Additional info

Anti-TdT antibody labels normal cortical thymocytes and primitive lymphocytes. Anti-TdT antibody detects an enzyme found in the nucleus of normal hematopoietic cells, normal cortical thymocytes and in the cytoplasm of megakaryocytes of the bone marrow. TdT expression is seen in over 90% of acute lymphoblastic lymphoma/ leukemia cases with the exception of pre-B-Cell ALL. TdT expression is not seen in normal mature T-or B-lymphocytes. Anti-TdT is positive for approximately one third of all cases of chronic myeloid leukemia, making it a good indicator of better response to chemotherapy.

Motea EA, et al. Biochimica et Biophysica Acta. 2010; 1804:1151-66

- Stauchen JA, et al. Int J Surg Pathol. 2003; 11:21-4
- 3. Suzumiya J, et al. J Pathol. 1997; 182:86-91
- 4. Arber DA, et al. Am J Clin Pathol. 1996; 106:462-8

5.

1.

2

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com